## **COLUMN USER GUIDE**

# for Agilent Normal-Phase and HILIC Columns

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This booklet provides general information for all Agilent ZORBAX, Agilent Pursuit, Agilent Polaris, and Agilent InfinityLab Poroshell normal-phase and HILIC columns. Initial sections include information for both normal-phase and HILIC columns. Because some method parameters and column care instructions are different for HILIC versus normal-phase, those sections are broken out separately in the last few pages. For additional information about your specific column or family, see **[agilent.com/chem/columnchoices](https://www.agilent.com/chem/columnchoices)**

## Getting Started

A QC Column Performance Report, including a test chromatogram, comes with every Agilent column. The QC test instrument is modified (optimized) from a standard instrument to minimize extra column volume, so results may vary from the system used in your lab. This allows a better evaluation of the column and assures a more consistent product.

## Installation

- The direction of flow is marked on the column.
- **1.8 µm columns (ZORBAX RRHT, ZORBAX RRHD) can only be operated in the flow direction marked on the column.**
- Agilent recommends using InfinityLab Quick Connect fittings (5067 5961 for 0.075 × 105 mm assembly) or Quick Turn Fittings (5067-5966) for high pressure column connections. Polyketone fittings (5042-8957) may also be used for pressures up to 600 bar. Be sure to order the right capillary connection for your fitting. Learn more at [www.agilent.com/chem/infinitylabfittings](http://www.agilent.com/chem/infinitylabfittings)
- Additional care must be taken to avoid damage to PEEK or PEEK-lined columns during installation. Combined compression and rotation may cause damage, therefore fittings without ferrules such as PEEK Finger-tight Fittings are not recommended. Fittings should not be overtightened, and should be swaged to the correct depth to avoid damage.

## Important Safety Considerations for all LC Columns

- All connection points in LC systems are potential sources of leaks. Users should be aware of the toxicity or flammability of their mobile phases.
- Because of the small particle size, dry column packings are respirable. Columns should only be opened in a well-ventilated area.
- Adhere to operating pressure limits noted for each column (see chart on next page). Exceeding these limits will compromise chromatographic performance and could be unsafe.

## Maximum Operating Pressures – Columns up to 9.4 mm id

Column type	Particle size	Pressure limit
Agilent ZORBAX NH2	5 µm, 7 µm	400 bar (6,000 psi)
Agilent ZORBAX HILIC Plus	3.5 µm	
Agilent ZORBAX Rx-SIL	5 µm	
Agilent ZORBAX SIL	5 µm, 7 µm	
Agilent ZORBAX CN	3.5 µm, 5 µm,	
Agilent ZORBAX Eclipse XDB-CN	7 µm	
Agilent ZORBAX Rapid Resolution High Definition (RRHD) HILIC Plus	1.8 µm	1,200 bar (17,000 psi)
Agilent ZORBAX RRHT Rx-Sil	1.8 µm	600 bar (9,000 psi)
Agilent Polaris NH2	3 µm, 5 µm	400 bar (6,000 psi)
Agilent Polaris Si-A	3 µm, 5 µm, 10 µm	400 bar (6,000 psi)
Agilent Pursuit XRs-Si	3 µm, 5 µm, 10 µm	400 bar (6,000 psi)
Agilent InfinityLab Poroshell 120 EC-CN	2.7 µm	600 bar (9,000 psi)
Agilent InfinityLab Poroshell 120 HILIC		600 bar (9,000 psi)
Agilent InfinityLab Poroshell 120 HILIC-Z		600 bar (9,000 psi)
Agilent InfinityLab Poroshell 120 HILIC-OH5		400 bar (6,000 psi)

### Column Usage Tips

- Modern columns are robust and are designed to operate for long periods under normal chromatographic conditions. You can maximize column lifetime by running it within specifications. Always review the specifications before putting in place a final method.
- While generally not harmful to the column, reverse flow should be avoided.
- Always use high purity reagents and chromatography grade solvents or better to prepare your mobile phase.
- Disassembling a column will degrade column performance.
- New columns contain a mixture of organic solvents. See your QC Performance Report for the solvent composition in your column. Initially, care should be taken not to pass any mobile phase through the column that may cause a precipitate or an emulsion to form.

- Agilent normal-phase columns are compatible with water and all common organic solvents.
- The use of a guard column is recommended to protect your column and increase its lifetime.
- Columns should not be maintained at high pH or elevated temperature when not in use.
- Avoid use of this column outside recommended pH range for column phase (see next page). Expect reduced lifetime when operating outside the recommended pH and temperature ranges.

### Tips for Getting the Best Chromatographic Results

- Optimize your instrument by minimizing tubing lengths between components to reduce extra column volume and band broadening. Use 0.12 mm id red tubing for normal operations. For narrow bore (2.1 mm or less) high resolution columns, such as InfinityLab Poroshell 120 or ZORBAX RRHD, use 0.075 mm id black tubing on UHPLC systems. Learn about capillary options at [agilent.com/chem/lccapillaries](https://www.agilent.com/chem/lccapillaries).
- When using normal-phase columns, ensure that your LC instrument's rotor seals, pump seals, tubing, and all other components are compatible with normal-phase solvents.
- Make sure your data collection rate is optimized for your column. Use a higher collection rate for Fast LC columns (InfinityLab Poroshell 120, RRHT, and RRHD).
- Use Agilent certified lamps in your LC instruments for best performance.
- Use sample filtration or other sample prep as appropriate for your sample. Learn more at [agilent.com/chem/sampleprep](https://www.agilent.com/chem/sampleprep).

## Working With Normal-Phase Columns

Every normal-phase column is tested before shipment and shipped in a normal-phase test eluent, with the exception of ZORBAX Eclipse XDB-CN, InfinityLab Poroshell 120 EC-CN, InfinityLab Poroshell 120 HILIC, InfinityLab Poroshell HILIC-OH5, InfinityLab Poroshell HILIC-Z, and ZORBAX HILIC Plus columns which are shipped in reversed-phase solvents. Care should be taken to make sure the column has been properly equilibrated prior to use. This will ensure reproducibility from analysis to analysis and help prevent retention time drifting.

### Normal-Phase Column Operating Parameters: pH and Temperature

Phase	Recommended pH range	Maximum operating temperature	Suggested operating temperature
Agilent ZORBAX NH2	pH 2.0 to 7.5	60 °C	30 °C
Agilent ZORBAX Rx-SIL Agilent ZORBAX SIL Agilent InfinityLab Poroshell 120 HILIC*† Agilent ZORBAX HILIC Plus**†	pH 1.0 to 8.0	Limited only by the temperature limits of the mobile phase	
Agilent ZORBAX CN	pH 2.0 to 8.0	60 °C	30 °C
Agilent ZORBAX Eclipse XDB-CN*	pH 2.0 to 9.0	60 °C	30 °C
Agilent Polaris NH2	pH 2.0 to 8.0	60 °C	30 °C
Agilent Polaris Si-A	pH 1.5 to 10.0	60 °C	30 °C
Agilent Pursuit XRs Si	pH 2.0 to 8.0	60 °C	30 °C
Agilent InfinityLab Poroshell 120 EC-CN*	pH 2.0 to 8.0	60 °C	30 °C

\* Shipped in reversed-phase solvents.

† See the *Working with HILIC* section for more information.

Operating at extreme ends of pH and temperature ranges will have a significant impact on column lifetime.

**Note:** All silica-based packings have some solubility in pH >6 aqueous mobile phases. When using silica-based columns at pH >6, best column lifetime is obtained at lower temperatures (40 °C Max) using low buffer concentrations in the range of 0.01 to 0.02 M.

### Normal-Phase Column Method Development Notes

The stationary phase is polar in nature and is best used with nonpolar mobile phases such as methylene chloride/hexane or isopropanol/hexane mixtures. Increasing the amount of the polar component in these mixtures typically reduces the retention time of the sample.

Column type	Method development notes
Agilent ZORBAX NH2 Agilent Polaris NH2	Retention decreases as the polarity of the mobile phase increases.  Many column volumes of mobile phase are necessary for equilibration to occur after changes in mobile phase composition. The number of column volumes required depends on the system in use and can be as high as 50 to 100 column volumes.
Agilent ZORBAX Rx-SIL Agilent ZORBAX SIL Agilent Pursuit XRs Si Agilent Polaris Si-A	Unbonded silica columns are very sensitive to trace amounts of water in solvents and sample, so care must be taken to maintain consistency for method reproducibility. When gradient elution is employed, both primary and secondary solvents should be modified with alcohol or acetonitrile and at least 30 columns of solvent should be allowed to flow through the column after completion of each run.  When switching between solvents, first rinse the column with a mutually miscible solvent such as isopropyl alcohol.
Agilent InfinityLab Poroshell 120 EC-CN Agilent ZORBAX Eclipse XDB-CN Agilent ZORBAX CN	Agilent InfinityLab Poroshell 120 EC-CN and ZORBAX Eclipse XDB-CN columns are shipped with reversed-phase solvents, but can also be used in normal phase applications after flushing with isopropanol.  Retention decreases as the polarity of the mobile phase increases.  The rapid equilibration of this bonded phase with mobile phase often makes it more convenient to use this column than a silica column, especially for gradient elution.

## **Cleaning Your Column/Extending Column Life: Normal-Phase Columns**

For normal-phase, we recommend use of organic solvents. Use at least 20 column volumes of each solvent, which is 50 mL if your column is the traditional analytical 4.6 × 250 mm. Try these solvents in order of increasing strength:

1. 50% methanol: 50% chloroform
2. 100% ethyl acetate

Cleaning a column used in normal-phase mode may depend on the sample type. If these solvents do not work, contact Agilent so we can recommend a solvent that may be more effective with your application and sample matrix.

### **Storage Recommendations: Normal-Phase Columns**

Normal-phase columns can typically be stored in their own mobile phase for long periods. Before storing the column, the end-fittings should be tightly capped with end-plugs to prevent the packing from drying out.

## **Working With HILIC Columns**

HILIC is a method used for the separation of polar analytes inadequately retained on typical reversed-phase columns. Agilent ZORBAX Rapid Resolution High Definition (RRHD) HILIC Plus columns and InfinityLab Poroshell 120 HILIC columns are non-bonded silica columns, shipped in acetonitrile:water or acetonitrile and ready to use for HILIC separations. InfinityLab Poroshell 120 HILIC-OH5 columns contain a polyhydroxy phase bonded onto superficially porous silica particles, and are shipped in acetonitrile. InfinityLab Poroshell 120 HILIC-Z columns contain a zwitterionic phase bonded onto superficially porous silica particles, and are shipped in acetonitrile:water. Like normal-phase columns, HILIC columns require more equilibration than reversed-phase columns. Note that ZORBAX Rx-Sil, ZORBAX Sil, and ZORBAX NH2 columns can also be used for HILIC applications. To change these columns to HILIC mode, flush first with isopropyl alcohol.

## HILIC Column Operating Parameters

Column name	pH Range	Maximum temperature	Suggested temperature
Agilent ZORBAX Rapid Resolution High Definition (RRHD) HILIC Plus Agilent ZORBAX HILIC Plus HPLC Agilent InfinityLab Poroshell 120 HILIC	0 to 8.0	60 °C	30 °C
Agilent InfinityLab Poroshell 120 HILIC-OH5	1.0 to 8.0	45 °C	30 °C
Agilent InfinityLab Poroshell 120 HILIC-Z	2.0 to 12	80 °C	30 °C

Operating at extreme ends of pH and temperature ranges will have a significant impact on column lifetime.

### Operating Guidelines: HILIC

- Agilent HILIC columns are shipped containing acetonitrile:water or acetonitrile.
- Care should be taken not to exceed column pressure limits, especially when changing mobile phase compositions.
- Agilent HILIC columns are compatible with water and all common organic solvents.
- Buffer should be mixed in small amounts, and refreshed often to prevent microbial growth.
- Eluents containing volatile acids (formic, acetic) or bases (ammonia) should be sealed tightly.
- All silica has some solubility in pH >6 aqueous mobile phases. Solubility of the silica is also increased at elevated temperatures. Using the ZORBAX RRHD HILIC Plus or InfinityLab Poroshell 120 HILIC columns above pH 6 and 40 °C will reduce the column lifetime.

## HILIC Method Development Notes

Column	Notes
Agilent ZORBAX Rapid Resolution High Definition (RRHD)	Most common HILIC applications are for the analysis of polar compounds that are unretained or poorly retained on reversed-phase columns.
Agilent HILIC Plus	It is best to equilibrate the column with 30-40% water in acetonitrile before use.
Agilent InfinityLab Poroshell 120 HILIC, HILIC-Z, and HILIC-OH5	Equilibration may require 20-50 column volumes. A typical mobile phase is acetonitrile:water with an acetate or formate buffer at 5-10 mM. Increasing buffer concentration can improve peak shape and retention. Decrease the aqueous/buffer component to increase retention.
Agilent InfinityLab Poroshell 120 HILIC-Z, and HILIC-OH5	<b>Suggested initial screening conditions:</b> A: 10 mM ammonium acetate B: (Acetonitrile:100 mM ammonium acetate, 90:10)* or A: 10 mM ammonium formate B: (Acetonitrile:100 mM ammonium formate, 90:10)* * For 10 mmol final concentration Adjust mobile phase pH as necessary Gradient: 90% B to $\geq 60\%$ B over 10 minutes

## **Cleaning Your Column/Extending Column Life: HILIC Columns**

Samples should be filtered before injecting on the column. This is particularly important for the ZORBAX RRHD HILIC Plus column, as the column inlet frit is nominally 0.5  $\mu\text{m}$  and samples should be filtered through a 0.2  $\mu\text{m}$  sample filter. If solvent flow appears to be restricted (unusually high column backpressure), check first to see that solvent flow is unobstructed up to the column inlet. If the restriction is prior to the column, replace the appropriate piece of tubing or filter that is plugged. If the column is plugged, InfinityLab Poroshell 120 HILIC can be backflushed with 20-30 column volumes of mobile phase. The RRHD HILIC column cannot be backflushed. To remove strongly-retained materials from the column, flush the column with a strong solvent which in HILIC mode, is water.

InfinityLab Poroshell HILIC-OH5 columns may be regenerated by flushing with 50/50 acetonitrile:50 mM ammonium acetate at a low flow rate for at least 3 hours. Then equilibrate with mobile phase.

InfinityLab Poroshell HILIC-Z columns may be regenerated by flushing with 50/50 acetonitrile:10 mM ammonium acetate at 20% of normal operating flowrates for at least 3 hours. Then equilibrate with mobile phase.

### **Storage Recommendations: HILIC Columns**

Acetonitrile:water (90:10) is recommended as the long-term storage solvent for HILIC columns. Before storing the column, tightly cap the end fittings with the end plugs to prevent the packing from drying out. Columns may be safely stored for short periods in most HILIC mobile phases. However, to protect equipment, it is best to remove salts from the instrument and column by purging the column with the same mobile phase without the buffer (for example, using 90:10 ACN:H<sub>2</sub>O to remove a 90:10 ACN:0.01M formate buffered mobile phase). Re-equilibration is faster with the original mobile phase when using this approach, but several (3 to 6) injections should be made to verify column equilibration.

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# LC Column and Sample Prep NAVIGATOR

## Find the Right Column for Your Separation

*Try the LC Columns Navigator*

With just a few mouse clicks, you can:

- Explore columns based on method parameters, compound, and USP method
- Find sample prep guidance and method optimization tips
- Get selection assistance from Agilent chromatography experts
- Obtain technical support contact information

Best of all, you can access the LC Columns Navigator anywhere – in your lab or in the field – from your smart phone or tablet

Get started now at [agilent.com/chem/navigator](https://www.agilent.com/chem/navigator)

## Agilent Ordering Information

For more information on our products and services,  
visit our web site at **agilent.com**

For technical support and local information,  
visit **agilent.com/chem/columnsupport**

To place an order,  
visit **agilent.com/chem/wheretobuy**

Agilent offers a complete line of sample preparation products to  
support LC and LC/MS applications.

The Agilent Bond Elut SPE and Captiva Filtration Sample Prep  
family of products offer the widest range of solutions for every  
level of sample cleanliness to help you increase throughput and  
enhance the quality of your data.

Learn more at **agilent.com/chem/sampleprep**

This information is subject to change without notice.

Agilent Technologies, Inc. 2014, 2018

Printed in Canada, January 24, 2018

820000-998